
# X6A end-station in 2002 (first sample)



- 1. Detector: CCD ADSC Q210
- 2. Diffractometer : (worm drive)  $\Theta$  spindle axis Crystal Logic with kappa phi orienter.
- 3. 45 degrees microscope with BNC CCD and monitor.
- 4. Goniometer head manually adjusted



### X6A end-station in 2005



- 1. ALS auto-mounter with 64 samples dewar
- 2. motorized XYZ (mini kappa)
- 3. "in-line" visualization with motorized zoom and focus + video server Using Actuated prism with pin hole.




#### X6A end-station in 2009



- 1. Detector: CCD ADSC Q270
- 2. single axis air bearing omega spindle
- 3. "in-line" visualization with motorized zoom and focus with fixed mirror (10 deg)
- 4. Large aperture collimating slits







<u>Detector</u>	Q210 (Thomson)	Q270 (Kodak)	
	210 mm	270 mm	<u>area * 1.65</u>
CCD QE(550 nm)	<0.25	<b>→</b> 0.7	efficiency * 3
Full well cap.	270 K	550 K	<u>dynamic range * 2</u>
Readout (full)	1 s	1 s	
PSF FWHM	90 um	90 um	

Rotation axis Worm drive Air bearing 6 deg/sec 180 deg/sec

>>>>> centering > 12times faster (increased screening speed by 2)

Collimating slits 1 mm 2 mm

Easy location of full beam, facilitate table / mirror alignment
(wish scannable slits with real time Io measurement)

# Helium flow compatibility

Slits/Ionization Chambers bar (60 + 40 cm): increase flux at sample at LE (\*4)

### Near in-line sample visualization


Robust, no motions, but motors. High magnification (FOV~ 240 wide)

## X6A end-station in 2013 (10 000 samples later)



Robust auto-mounter
Robust centering
Robust lid
Remote friendly
Well equipped (vortex,
AVT, annealing, ...)

Possible improvements (that would be used at NSLS II) K-B mirror to increase flux density by > 10 High capacity auto-mounter dewar holding 7-10 pucks UHR GigE camera, 14 MP, motionless zoom 1 to 5







X6A: current and past staff members

NSLS: staff (Sci. Tech. Mech. IT.)

